



CIVCO
MEDICAL SOLUTIONS



Safe. Efficient. Guiding Patient Care.

ULTRA-PRO ™

1

attach bracket

Attach custom reusable non-sterile bracket to transducer by aligning locating features. CIVCO manufactures brackets for Aloka, Biosound Esaote, BK Medical, Fukuda, GE Healthcare, Hitachi, Kontron, Medison, Philips, Shimadzu, Siemens, SonoSite, Terason, Toshiba, Ultrasonix and ZONARE.

2

safe. efficient. go

Six steps for clinical confi

4

gauge selection

Easy-to-read gauge sizes on the inserts makes it simple to identify and alter gauge sizes in a darkened ultrasound suite. The guide also features a large funnel for instrument insertion and will accept a full range of sizes including: 8.5FR, 14-23GA (19GA not available).

5

cover up

CIVCO makes it easy to reduce infection rates and cross-contamination with a telescopically-folded CIV-Flex™ cover. Expand your sterile field to 58" by covering the transducer and cord with the sterile latex-free cover.



needle guide assembly

Use proper sterile technique to attach needle guide onto bracket and lock into place with a single click.

Guiding patient care.

vidence in biopsy procedures.

place guide insert

After selecting the appropriate gauge size, the insert slides easily into the needle guide and is ready for use providing a fixed path that directs instruments according to on-screen biopsy guidelines.



efficient quick-release

The Ultra-Pro II features a large tab for improved quick-release function, allowing easy detachment of the needle from the transducer assembly.

Safe. Efficient. Guiding Patient Care.

Ref. #	Product Description	Qty.
--------	---------------------	------

Needle Guide Kits - Sterile procedure kits.

610-699	Sterile Ultra-Pro II™	24
610-579	Sterile Ultra-Pro II with 14 x 61cm (5.5" x 24") CIV-Flex™ cover (3D)	24
610-608	Sterile Ultra-Pro II with 14 x 91.5cm (5.5" x 36") telescopically-folded CIV-Flex cover (3D)	24
610-904	Sterile Ultra-Pro II with 14 x 147cm (5.5" x 58") telescopically-folded CIV-Flex cover (3D)	24
610-609	Sterile Ultra-Pro II with 17.8 x 91.5cm (7" x 36") telescopically-folded CIV-Flex cover (3D)	24
610-917	Sterile Ultra-Pro II with 17.8 x 147cm (7" x 58") telescopically-folded CIV-Flex cover (3D)	24

Transducer Covers - Protective sheath used to prevent cross-contamination.

610-1000	Sterile 10.2 x 147cm (4" x 58") CIV-Flex cover	24
610-542	Sterile 14 x 91.5cm (5.5" x 36") CIV-Flex cover (3D)	24
610-575	Sterile 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)	24
610-1010	Sterile 17.8 x 147cm (7" x 58") CIV-Flex cover (3D)	24

System Drapes - Protective sheath used to protect equipment.

610-1037	Sterile 20 x 28cm (7.9" x 11") polyethylene system drape	24
610-025	Sterile 76 x 76cm (30" x 30") polyethylene system drape	24

CIVCO manufactures brackets for Aloka, Biosound Esaote, BK Medical, Fukuda, GE Healthcare, Hitachi, Kontron, Medison, Philips, Shimadzu, Siemens, SonoSite, Terason, Toshiba, Ultrasonix and ZONARE.

For bracket availability call 319.656.4447 or visit www.civco.com.

Safe. Efficient. Guiding Patient Care.

Introducing the Ultra-Pro II™, a needle guide that incorporates more flexibility and advanced features to increase productivity during ultrasound-guided procedures. The single-use design of the guide reduces the risk of cross-contamination and increases throughput. The guide also offers physicians reduced technique variability, providing a shorter learning curve and reduced procedure time.



Accurate Placement

Consisting of a two-part system of a reusable bracket and a disposable snap-on needle guide, the Ultra-Pro II provides a safe, reliable and efficient tool for ultrasound-guided procedures including catheter placement, core tissue biopsy, drainage procedures, fine needle aspirations and regional anesthesia. The Ultra-Pro II holds the needle as it is advanced, providing a smooth insertion, while allowing you to see what you puncture.



Prevention
Precision
Positioning

- **Prevention:** Transducer and equipment covers protect against cross-contamination providing patient and staff safety.
- **Precision:** Needle guides are designed to reduce technique variability in order to rapidly deliver accuracy and stability for puncture procedures.
- **Positioning:** Needle guides deliver accurate positioning during puncture, providing confidence in hitting the target.

Safe. Efficient. Guiding Patient Care.

ULTRA-PRO ™

Placing an Order

Contact a CIVCO customer service specialist to place an order or request a product sample.
For international distribution, visit CIVCO Worldwide at WWW.CIVCO.COM

Toll-free within the U.S. and Canada: 800.445.6741

Toll-free Facsimile: 877.329.2482

Direct Dial: 319.656.4447

Facsimile: 319.656.4451

E-mail: info@civco.com

WWW.CIVCO.COM

